metrology is our business, accuracy is our standard





Technical Specification



www.amtest-tm.ro

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.ro 1

Departament Vanzari: +40 73 071 4926 Iudmila.valcea@amtest-tm.ro



Content

Content
1. Equipment Appearance:
2.Equipment Information:
2.1 Application:
3.Technical data:
3.1Basic information:
3.2 Temperature:
3.3 Humidity:
3.4 Sample limitation and test method:
4.Equipment structure:
5.Temperature and Humidity System:
6.Refrigerant System:
7.Control System:
8.Main Sapre Parts:
9.Safty Devices:10
10.Ambient Environment:

www.amtest-tm.ro

+40 21 414 0401 office@amtest-tm.ro **Departament Vanzari:** +40 73 071 4926 Iudmila.valcea@amtest-tm.rc



1. Equipment Appearance (For reference)

2. Product description

2.1 Product applications

The constant temperature and humidity test chamber is widely used in scientific research, research institutes, quality inspection institutes and other industries to do warm and humid heat tests of electronic and electrical products, materials, parts, equipment, etc., alternating heat test and constant temperature test, etc. Also could do high and low temperature routine tests, low temperature storage, in order to evaluate the performance of the sample under given environmental conditions.

2.2 Product Positioning

Apply to aerospace, aviation, electronics, automotive, battery and other products and quality inspection institutes, research institutes, universities and other test units, providing virtual space to simulate the real environment, verify product inspection, research and development results, test chamber is a powerful assistant to shorten the development cycle and improve product quality and reliability.

3. Technical indicators

3.1 Basic information

Model	SMC-225-CC
Temperature control range	-40°C~150°C
Humidity control range	20%RH~98%RH
Condenser	Air-cooling
Refrigerant	R449A
Interior size(mm)W*H*D	600x750x500
Outer size(mm) W*H*D	800x1870x1315
Volume(L)	225L
Operating ambient temperature	$+5^{\circ}C \sim 35^{\circ}C$
Power supply	220V AC, 50Hz single phase
Controller	Sanwood programmable controller with USB port

3

www.amtest-tm.ro

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.ro **Departament Vanzari:** +40 73 071 4926 ludmila.valcea@amtest-tm.ro



3.2 Temperature test

Temperature range	-40°C∼150°C
Temperature fluctuation	≤±0.5°C
Temperature deviation	l≤±2°C
Temperature uniformity	≤1.5°C
Temperature resolution	0.01°C
Heating rate	$25^{\circ}C \rightarrow +100^{\circ}C$ /within 25 mins (with standard load)
Cooling rate	150.0°C~25.0°C Cooling rate 2.0~3.0°C/min
	25.0°C~-40.0°C Cooling rate 1.0~2.0°C/min
Standard load	20kg aluminum sheet, 200W heat load



3.3 Humidity test

Humidity range	20%R.H \sim 98%R.H when temperature from 20°C ~ 85°C
	(in above map of blue area)
Humidity deviation	$\leq \pm 2\%$ R.H
Humidity fluctuation	±2%R.H
Humidity resolution	0.1%R.H
Verify conditions	1) Condition 1: -40°C±2 °C
	2) Condition 2: 22°C±2 °C, 95%±2%RH
	3) Condition 3: 105°C±2 °C
	4) Condition 4: 22°C±2 °C, 95%±2%RH
	5) Change over time from condition 1 to 2, condition 2 to

4

www.amtest-tm.ro

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.ro **Departament Vanzari:** +40 73 071 4926 Iudmila.valcea@amtest-tm.ro



condition 3, condition 3 to 4 within 60 ± 5 minutes.

3.4 Sample limit and test method

Prohibitions	Testing and storage of explosive, flammable, volatile materials
	Testing and storage of corrosive substances
	Testing or storage of biological samples
	Test and storage of strong electromagnetic emission source samples
	Testing and storage of radioactive material samples
	Testing and storage of samples of highly toxic substances
	Testing and storage of samples that may produce highly toxic substances
	during testing or storage
Test	GB/2423.1-2008 (IEC60068-2-1:2007) low temperature test method AB.
standard	GB/T5170.5-2008 damp heat test equipment.
	GJB150.4 (MIL-STD-810D) low temperature test method.
	GB2423.3-93 (IEC68-2-3) Test Ca: Constant damp heat test method.
	GB2423.4-93 (IEC68-2-30) Test Db: Alternating Damp Heat Test Method

4. Machine structure

	-
Structure	One piece assembled type
Inner chamber	SUS#304 heat-resistant and cold-resistant stainless steel plate
material	(1.2mm); full seamless welding
Internal structure	SUS304 (2mm) stainless steel reinforcement
strengthening	
Outer chamber	Electrolytic steel sheet 1.2mm, pickling phosphating high-grade
material	powder baking varnish
Insulation material	Germany Bayer refractory grade high strength PU polyurethane
	foam insulation material 100mm and ultra-fine glass fiber 10mm
Door edge	Double-layer high-tension silicone rubber seal, temperature resistant
	-90~180°C, lifespan up to 15 years
Observation	multi-layer hollow tempered glass belt with automatic defrosting
window	function, which can guarantee the frost-free and condensation of the
	glass surface during any test. Internal lighting over the window.
Sample rack	Stainless steel sample holder 2 units, height is adjustable
	load-bearing (uniform): 20kg/unit
Moving and	4 high load-bearing pulleys and PU horizontal angle wheels at the
position mode	bottom for moving and fixing the equipment
Cable port	Two cable ports, on the chamber's left and right side,φ100mm.
	with screwed plastic cover, silicone plug
Floor bearing	$\leq 100 \text{kg/m2}$ (uniform load)

www.amtest-tm.ro

5

Office: Bucuresti, Romania Cara Anghel 9. bl. C56. ap.8

+40 21 414 0401 office@amtest-tm.rc **Departament Vanzari:** +40 73 071 4926 udmila.valcea@amtest-tm.ro



Circulating motorStainless steel extended shaft circulating motor ensures long-term operation and sufficient air volume operationCirculating windThe multi-wing centrifugal circulating wind wheel is used to strengthen the shaft and aluminum alloy to make high and low temperature resistant rotating blades, so as to achieve forced convection and effectively avoid looping dead angles.Circulating air ductThe temperature-adjusting and conditioned air duct is designed as a double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		
operation and sufficient air volume operationCirculating windThe multi-wing centrifugal circulating wind wheel is used towheelstrengthen the shaft and aluminum alloy to make high and lowtemperature resistant rotating blades, so as to achieve forcedconvection and effectively avoid looping dead angles.Circulating air ductThe temperature-adjusting and conditioned air duct is designed as adouble air duct, which is connected to the studio but isolated. Thewind path from up to down and returning. The partition plate isformed by cold-bending processing of high-quality stainless steelplates, and adjustable louvers are used at the air outlet. Indirectheater, saturated humid air inlet, refrigeration dehumidificationevaporator and circulating blast wind wheel are arranged in thetemperature regulation air passage.	Circulating motor	Stainless steel extended shaft circulating motor ensures long-term
Circulating wind wheelThe multi-wing centrifugal circulating wind wheel is used to strengthen the shaft and aluminum alloy to make high and low temperature resistant rotating blades, so as to achieve forced convection and effectively avoid looping dead angles.Circulating air ductThe temperature-adjusting and conditioned air duct is designed as a double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		operation and sufficient air volume operation
wheelstrengthen the shaft and aluminum alloy to make high and low temperature resistant rotating blades, so as to achieve forced convection and effectively avoid looping dead angles.Circulating air ductThe temperature-adjusting and conditioned air duct is designed as a double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.	Circulating wind	The multi-wing centrifugal circulating wind wheel is used to
temperature resistant rotating blades, so as to achieve forced convection and effectively avoid looping dead angles.Circulating air ductThe temperature-adjusting and conditioned air duct is designed as a double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.	wheel	strengthen the shaft and aluminum alloy to make high and low
convection and effectively avoid looping dead angles.Circulating air ductThe temperature-adjusting and conditioned air duct is designed as a double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		temperature resistant rotating blades, so as to achieve forced
Circulating air duct The temperature-adjusting and conditioned air duct is designed as a double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		convection and effectively avoid looping dead angles.
double air duct, which is connected to the studio but isolated. The wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.	Circulating air duct	The temperature-adjusting and conditioned air duct is designed as a
wind path from up to down and returning. The partition plate is formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		double air duct, which is connected to the studio but isolated. The
formed by cold-bending processing of high-quality stainless steel plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		wind path from up to down and returning. The partition plate is
plates, and adjustable louvers are used at the air outlet. Indirect heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		formed by cold-bending processing of high-quality stainless steel
heater, saturated humid air inlet, refrigeration dehumidification evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		plates, and adjustable louvers are used at the air outlet. Indirect
evaporator and circulating blast wind wheel are arranged in the temperature regulation air passage.		heater, saturated humid air inlet, refrigeration dehumidification
temperature regulation air passage.		evaporator and circulating blast wind wheel are arranged in the
		temperature regulation air passage.

5. System introduction

Heating wire	High-quality explosion-proof nickel-chromium alloy heating wire (high
	resistivity, small temperature coefficient of resistance, small
	deformation at high temperature and not easy to embrittlement, self-
	heating temperature up to 1000-1500°C, long service life) rapid heat
	exchange, no hysteresis
Heating wire	The solid state relay is used as a heating actuator, and there is no large
control	current fluctuation and impact phenomenon, and the operation is stable.
Heating wire	The heating wire is provided with anti-dry protection to prevent the
protection	heater from continuously burning after the circulation fan stops for some
	reason, causing the heater itself to burn out or other accidents.
Humidification	Steam humidification method: using electronic parallel mode micro-
mode	motion humidification system
Humidifying	All stainless steel embedded humidification tube with anti-dry explosion
heating pipe	protection protector
Humidification	Humidification is rapid, saves water, saves electricity, and allows test
system	products to heat up. The humidification and dehumidification system are
advantage	completely independent, no need for extra drainage, faster than
	traditional surface humidification (water tray), high control precision, no
	scale pollution such as scale and scale, good low humidity performance,
	water level observation window, and easy cleaning
water for	Distilled water, pure water or deionized water (resistivity greater than
humidity	500 $\Omega \cdot m$) to be used for humidity system (provided by the user)
Water storage	Drawer type water tank 1 unit, 15L, located in front of the machine,

www.amtest-tm.ro

6

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.ro **Departament Vanzari:** +40 73 071 4926 udmila.valcea@amtest-tm.rc



device	under the door.

6. Refrigeration System

0	
Compressor	France Tecumseh fully enclosed compressor
Evaporator	High-efficiency components adopts a slope type evaporator (AC&R
	compound spoiler aluminum fins)
Condenser	Air-cooled system for equipment easy movement, etc.
Heat exchanger	SWEP plate type refrigerant cold and heat exchange design, making
	higher efficiency compared with traditional internal spiral
Energy-saving	Adopting throttle electronic expansion valve
device	1. The active control of the refrigeration system is realized, and the
	fixed proportional adjustment of the original thermal expansion valve
	is not controllable. The output can be adjusted in advance and
	optimized for different modes and operating conditions.
	2. Due to the cyclic control feedback of the electronic expansion valve,
	the front end is a temperature-plus-pressure dual-sensor high-response
	direct control, which can provide the best evaporator liquid supply, so
	that the refrigeration system can achieve excellent cooling capacity in a
	wider working range. Output.
	3. Energy saving: the full range of electronic expansion valve self-
	adjustment + active adaptation to adjust the cooling capacity output,
	making the system more energy efficient.
	\$ 67
	4. Energy-saving design: Adopting PID + PWM principle of VRF
	(refrigerant flow control) technology to achieve low-temperature
	energy-saving operation (electronic expansion valve according to
	thermal energy conditions refrigerant flow servo control technology)
	low temperature working state, the heater does not participate in the
	work, through PID + The PWM regulates the flow rate and flow

www.amtest-tm.ro

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.ro 7

Departament Vanzari: +40 73 071 4926 Iudmila.valcea@amtest-tm.ro



	direction of the refrigerant, and regulates the three-way flow of the
	refrigerant pipe, the cold bypass pipe, and the hot bypass pipe to
	achieve automatic constant temperature to the working chamber. This
	method can reduce energy consumption by 30% under low temperature
	conditions. The technology is based on Danish Danfoss' ETS series of
	electronic expansion valves, which can be used to smoothly adjust the
	cooling capacity for different cooling capacity requirements, that is, to
	achieve the compressor cooling capacity adjustment when different
	cooling rate requirements are met.
Refrigerant	Original American DuPont Environmental Refrigerant R449A

7. Control system

Controller	Sanwood program controller
Display	640x480 dot matrix, 5.7 inch TFT color LCD display
Running mode	Program mode, fixed value method
Setting mode	Chinese and English menu (free choice), touch screen input
Program capacity	Editable program
	Quantity: Max 120
	Steps: Max 100
	Number of cycles: Max 999
	Program can be linked (link program serial number can be
	selected)
Setting range	Temperature: Adjust according to the operating temperature range
	of the equipment (upper limit +5°C, lower limit -5°C)
	Humidity: $(0 \sim 100)$ % RH (temperature and humidity test
	equipment)
Display resolution	Temperature: 0.01°C
	Time: 0.01min
	Humidity: 0.1% RH (temperature and humidity test equipment)
Communication	Ethernet, RS485, RS232 interface;
function (standard)	choose one from above three options;
Control mode	BTC balance temperature control mode + DCC (Intelligent
	Cooling Control) + DEC (Smart Electrical Control) (Temperature
	Test Equipment)
	BTHC balance temperature regulation and humidity control mode
	+ DCC (smart cooling capacity control) + DEC (smart electrical
	control) (temperature and humidity test equipment)
Curve recording	With battery-protected RAM, it can save the set value, sampling
function	value and sampling time of the device; the maximum recording

www.amtest-tm.ro

8

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm. **Departament Vanzari:** +40 73 071 4926 Iudmila.valcea@amtest-tm.rc



	time is 350 days (when the sampling period is 1.5min)	
Subsidiary function	Fault alarm and cause, processing prompt function	
	Power failure protection	
	Upper and lower temperature protection	
	Calendar timing function (automatic start and automatic stop	
	operation)	
	Self-diagnostic function	

8. Main parts

Name	Brand	Note
Controller	Sanwood	SA <mark>\</mark> ₩OODΞπ
Compressors	France Tecumseh	🧕 Tecumseh
Oil separator	Emerson	EMERSON
Heat exchanger	SANWOOD customized	SA <mark>\</mark> ₩OOD Ξ π
Condenser	SANWOOD customized	SA <mark>\</mark> ₩OODΞπ
Evaporator	SANWOOD customized	SAℕWOOD三市
Drying filter	Denmark Danfoss	Danfoss
Non-return valve	Denmark Danfoss	Danfoss
Solenoid valve	Denmark Danfoss	Danfoss
Condensing pressure switch	Denmark Danfoss	Danfoss
Leakage protection switch	France Schneider	Schneider
AC relay	France Schneider	Schneider
Thermal relay	France Schneider	Schneider
phase rotation relay	Switzerland Carlo Gavazzi	_//∧LŮ.佳牙
solid-state relay	Switzerland Carlo Gavazzi	_//\L₽. 佳牙
Temperature humidity sensor	Taiwan Thermoway	
Circulating motor	Taiwan Yuzheng	
Over temperature protection	South Korea RAINBOW	

www.amtest-tm.ro

omce: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.ro 9

Departament Vanzari: +40 73 071 4926 ludmila.valcea@amtest-tm.rc



9. Equipment safety protection device

Refrigeration	Compressor over-pressure protection
system	Compressor motor overheat protection
	Over current protection of compressor motor
Humidifying	Humidification tube dry burning protection
system	Abnormal water supply protection
	Abnormal drainage protection
	Water shortage abnormal protection
Test chamber	Adjustable test product over temperature protection
protection	Test space temperature fuse protection
	Air conditioning channel limit over temperature protection
	Ultra-high temperature protection inside the controller
	Fan motor overheat protection
	Chamber inside and outside pressure balance protection
Other	Total power phase sequence and phase loss protection
protection	Whole machine leakage protection
	Load short circuit protection

10. Equipment safe use conditions

Site	Leveled ground, good ventilation, no flammable, explosive or
	aggressive gas and powder
	No strong electromagnetic radiation source
	With drainage near the machine (within 2 meters)
	Ground load-bearing capacity: not less than 800kg/m2
	Leave space around the machine for maintenance
	A : not less than 30cm
	B : not less than 50cm
	C : not less than 70cm
	D : not less than 90cm
Power supply	Voltage allowable fluctuation range: AC (1±10%) 220V or
	AC (1±10%) 380V
	Frequency allowable fluctuation range: $(1 \pm 1\%)$ 50Hz
	Protective earthing wire grounding resistance is less than 4Ω
	Users are required to configure a device with a considerable capacity of
	air or power switch at the installation site, and this switch must be used exclusively for this device.

www.amtest-tm.ro

10

Office: Bucuresti, Romania Cara Anghel 9. bl. C56. ap.8

+40 21 414 0401 office@amtest-tm.ro **Departament Vanzari:** +40 73 071 4926 ludmila.valcea@amtest-tm.rc



Requirements	When the equipment is not working, the ambient temperature should be
for storage	kept within $0 \sim +35^{\circ}$ C (without freezing)
environment	When the ambient temperature is lower than 0 °C, the water remaining
	in the equipment should be drained to prevent the water in the pipeline
	from freezing and damaging the pipeline.
Other	Opening the door of the test chamber during the test will cause
	temperature or humidity fluctuations inside the box
	If the door is opened several times during the test or the door is opened
	for a long time or the test sample emits moisture, the heat exchanger of
	the refrigeration system may freeze and may not work properly.

11. Quality assurance

Warranty: 3-year warranty included, calculated from the date of production. Extension warranty is negotiable.

After-sale service:

1) the equipment is free of charge during the warranty period, and the end user is in compliance with the conditions of custody, use and installation rules.

2) Due to the failure caused by the quality, design defects and/or core components of the chamber, AMTEST-TM România will provide the technical support and solutions until the chamber runs well. In-site service and consumable components are not included.

3) Natural disasters, abnormal power, improper use, damage caused by improper maintenance, etc are not in the warranty scope.

Documents and accessories:

Technical information: product certificate, instruction manual, warranty card and etc.; **Packaging and shipping methods**

Packing: Shipping container that meets the requirements of QB/BWD008-2001 **Mode of transport**: by Sea, rail or land.

Training: Your operator can operate the machine skillfully according to the manufacturer manual. Video will be provided if required.

Installation: One piece assembled machine.

Maintenance: AMTEST-TM România will do the maintenance surcharge according to the manufacturer guidance.

Romania +40 2 21 9, bl. C56, ap.8 office **Departament Vanzari:** +40 73 071 4926 ludmila.valcea@amtest-tm.





www.a

Office: Bucuresti, Romania Cara Anghel 9, bl. C56, ap.8

+40 21 414 0401 office@amtest-tm.rc **Departament Vanzari:** +40 73 071 4926 ludmila.valcea@amtest-tm.ro